

02 Safety Precautions


When using the massage chair, basic precautions should always be followed to reduce the risk of electric shock, burns, fire or injury to persons. Please read the following precautions before using the massage chair.

- * Make sure the power plug is securely inserted into a proper electrical outlet to reduce the risk of short-circuit and fire.
- * Always unplug the massage chair from the electrical outlet immediately after using and before cleaning.
- * **DO NOT** use the massage chair in areas where the humidity is high such as a sauna or steam room.
- * **DO NOT** use the massage chair outdoors and avoid areas with direct sunlight.
- * If the electricity is off, turn the power switch to the OFF position and unplug the massage chair from the electrical outlet to prevent any damages to the chair when the electricity comes back on.
- * Never plug/unplug the power cord plug socket with wet hands.
- * Close and continuous supervision is necessary when the massage chair is used by, on or near children or disabled persons.
- * No more than one person should use the massage chair at a time.
- * The massage chair is intended for use as described in the Instruction Manual. Do not use any attachments other than those recommended by the manufacturer.
- * **DO NOT** attempt to stand or sit on the backrest, armrest or footrest in order to prevent accident, injury or damage to the chair.
- * Avoid falling asleep while using the massage chair.
- * Always keep the massage chair clean, and never insert any objects into or around the gaps of the massage chair.
- * **DO NOT** move the massage chair by pulling on the power supply cord.
- * Make sure there are no children, pets or other obstacles behind or underneath the massage chair.
- * If you find any rips or tears in the fabric of the massage chair contact the service department and stop using the chair immediately.

-
- * This massage chair is designed for household use and should not be used in non-residential settings.
 - * Please consult your doctor if you are in doubt as to whether or not you should use this massage chair.
 - * Please consult a doctor before using the massage chair if you are pregnant, undergoing any medical treatment, using a medical electronic device such as a pacemaker, or are experiencing any of the following: Malignant tumor, heart disease, back pain, abnormal or curved back bone, osteoporosis or acute illness.
 - * If you experience pain while using the massage chair, stop using the chair immediately and consult your doctor.
 - * Never use a massage function on swollen or inflamed parts of the body.
 - * DO NOT use the massage chair for more than 30 minutes at a time.
 - * DO NOT use the massage chair where aerosol (spray) products are being used or where oxygen is being administered.
 - * This massage chair is a non-professional product designed to provide a comfortable massage and should not be substituted for appropriate medical treatment.
 - * Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair, and the like.
 - * Temperatures sufficiently high to cause burns may occur regardless of the control setting. Do not use on an invalid or on a sleeping or unconscious person. Check the skin in contact with the heated area of the appliance frequently to reduce the risk of blistering.
 - * Burns may result from improper use.
 - * Do not crush-avoid sharp folds.
 - * Keep children away from extended foot support (or other similar parts.)
 - * DO NOT WET-DO NOT USE PINS.
 - * NEVER REMOVE COVER.

03 Chair Components





air bags

0 air bags
hammer

air bags
hammer

0 air bags



















air bags

04 Chair installing



When installing the chair, please leave 31.5" of space behind the back-rest to allow free reclining.
Unplug from outlet when installing parts of the massage chair.

Before installing, please check all the parts in the following list.


1	 Massage chair part	2	 Backrest	3	 Side panel
4	 Backrest pad	5	 Back massage strength-control pad	6	 Pillow pad
7	 Controller holder	8	 Power supply wire	9	 Flat head rivet
10	 Split pin	11	 Phillips screwdriver	12	 φ6*20*2 Nylon spacer
13	 M4*10 Pan head screws×4	14	 φ4*9*1 Gasket×4	15	 M6*55 Hex screws×4
16	 φ6*18*2 Gasket×4	17	 φ6 Spring washer×4	18	 Allen driver

Step 1: Remove the chair parts out of the two cartons:

- 1) Massage chair part
- 3) Side panel

Note: The parts and the serial number of them in the below installation are corresponding to the parts list.

Step 2: Connect pipes and terminals between the side panels and seat part as demonstrated in PIC 1.





Step 3: Install side panels.

- 1) Insert the side panels along the hook D to the holder C (on the seat) at 90 degrees as demonstrated in PIC 2.
- 2) Adjust the side panel to make the holes aligned, lock the hex screws as demonstrated in PIC 3.

Note:

- 1) Please be careful to avoid scratching during installation.
- 2) The hex screws, gaskets and spring washers are packed in a small bag marked "1".
- 3) Please always bear in mind to fix the side panels to the chair without backrest, meanwhile, lifting up the seat pad for easier to tighten the screws. Once you fix one of the side panels, please don't tighten the two screws immediately until you make sure the two screws are perfectly fitting to each holes.


4) Please fix the back screw prior to the front screw while you assemble one of the side panels.



Step 4: Repeat step 2 and step 3 to install the other side panel.

Step 5:


- Use the diagonal pliers to cut off the cable tie which binds the stopper of the backrest and take off the packing material of the backrest actuator, where both of them are located underneath the seat pad.



Step 6: Install the backrest.


Lift the backrest, then lower the backrest into the seat brackets. Securely lock the flat head bolts (A) into both sides (B) of the seat as demonstrated in PIC 5.

- A. Flat Head Bolt(on the backrest)
- B. Holder(on the seat)




Step 7: Attach the connecting rod (A) and the U stay fork of the reclining actuator (B) with a flat head rivet (C), and secure with a split pin (D), then connect the air pipe and the audio wire in both left and right sides as demonstrated in PIC 6.

- A. Connecting rod (on the backrest)
- B. U stay fork of reclining actuator (on the seat)
- C. Flat head rivet (part 10 in the list)
- D. Split pin (part 11 in the list)




Step 8: Insert the wire plugs into the seat as demonstrated in PIC 7.

Note: Insert black to black and grey to grey, otherwise damage to the chair may happen.




PIC 7




Step 9: Attach the backrest pad and the pillow pad.

- 1) Attach the backrest pad and the pillow pad using the magic tape.
- 2) Attach the backrest pad with the zipper and magic tapes.



PIC 8

Step 10: Attach the outer coat of the massage pillow that is easy to take off and could be cleaned by water.




PIC 9

Step 11: Install the controller.


Adjust the positioning holes of controller to match with the ones in the side panel as demonstrated in PIC 11. Lock the screws which are with flat pad when holes C', D', E', F' are corresponding to hole C, D, E, F.

Note: 1. The A and B in PIC 11 are corresponding to PIC 10's A and B.

2. The screws and the gaskets used in this step are packed in a small bag marked "3".




PIC 10




PIC 11


Step 12: Connect the wire plugs.




PIC 12



Step 13: Put the instruction manual into the small pocket on the backrest pad. The small pocket's position is shown in PIC 13.




Step 14: Music player can be put inside the remote controller's holder as demonstrated in PIC 14.




Note: Please use the nylon spacer (put in small bag marked "4") for proper adjustment if the gaps between the side panels and the backrest are different.

05 Preparation For Use

- To set up the chair
 1. To avoid obstruction to the back and front of the chair, please ensure when the chair reclines, there will be excess space free from any obstacles.
 2. There should be at least 10" for the leg rest and 31.5" for the backrest.
 3. To avoid damage to the floor, it is advisable to place a carpet or cushion under the chair.



- To move the chair
 1. Lift up the leg rest.
 2. Recline the chair to approximately 45 degrees.
 3. Move the chair around the house with a slight pull and push.



Caution:

- Before moving the chair, adjust the leg rest back into its original position.
- Do not move the massage chair when there is someone seated on it.
- Place a mat underneath the massage chair before moving to protect the floor from scratches.