

NEOMASK Storage and Care Instructions:

SKU: SG-NEO-SINGLE

To retain product quality, NEOMASK should be stored in a cool, dry place when not in use.

The NEOMASK neoprene material is hand washable, remove the filters before cleaning. Do not clean NEOMASK filters, when soiled, replace them.

Product should be handled with care when removing the Activated Carbon Cloth.

Wash Instructions:

- Remove carbon/external filters
- Hand wash separately in cool water
- Use mild detergent
- Hang to dry

NEOMASK Storage and Care Instructions:

SKU: SG-NEO-Filters

Non-woven filters should not be cleaned or reused. Discard soiled filter and replace with a new filter.